



Fetuin

Catalog Number SA-F001

Lot no. 404.1A

Product Fetuin glycoprotein standard

Source Bovine


Amount 0.5 mg

Storage Keep at 4°C

Molecular Weight 48000

Shipping: The product should be shipped overnight on cold packs.

Description Fetuin is a heavily glycosylated protein which exists in a variety of glycoforms containing bi-, tri-, and tetra-antennary oligosaccharides with variable sialylation. The most abundant N-glycans on the glycoprotein are the triantennary species [Hayes et al (1992), Townsend et al (1989)]. Fetuin is commonly used as a glycoprotein control for carbohydrate release and labeling.



Warranties and liabilities

QA-Bio, LLC warrants that the above product conforms to the specifications described herein. Should the product fail for reasons other than through misuse QA-Bio, LLC will, at its option, replace free of charge or refund the purchase price. This warranty is exclusive and QA-Bio, LLC makes no other warrants, expressed or implied, including any implied conditions or warranties of merchantability or fitness for any particular purpose. QA-Bio, LLC shall not be liable for any incidental, consequential or contingent damages.

This product is intended for *in vitro* research only.

updated 5/22/2007

SA-F001 Fetuin Specifications